



Southeast Scientific
770.426.6678

Customer	CHOA - Center for Advanced Pediatrics		
Address	2174 N Druid Hills Road NE		
City, State, Zip	Atlanta, GA, 30329		
Ph	404-852-0765	Contact	Mimi Le
	1 st PM <input checked="" type="checkbox"/>		2 nd PM <input type="checkbox"/>

Unit Preventative Maintenance – Refrigerator / Freezer

Date of Service: 9/19/24 Next Service Due: 3/19/25 Room # 5126
 Model # HCUCFS-0420 Serial # VWR-041805280116-1808 Mfgr VWR
 SES # 9424 Customer # OKOYE Tech CD
 NIST Meter Serial # T-281182 NIST Meter Calibration Due 6/26/26

	Pass	Fail	N/A	NEW		Yes	No	N/A
Door Seal	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Evap Drain Clear	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Hinges	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		CFM	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Battery	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Clean Condenser	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Gaskets	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Insulation Lines	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Handles	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					

Temperature									
Setpoint	<u>D 7</u>	Display	<u>NA</u>	NIST	<u>-26.6C</u>	Chart Recorder	<u>NA</u>	Adjust	<u>NA</u>

Compressor		
<u>115</u> Voltage	<u>0.5</u> Amperage	<input type="checkbox"/> Inverter Normal <input checked="" type="checkbox"/> No Inverter <input type="checkbox"/> Inverter Error (See Notes)

1.0 Unit Total Amperage NA High Alarm Set NA Battery Charge Circuit
0.9 Amps Rated NA Low Alarm Set GOOD Wiring Connections
115 Volts Rated

Technician's Notes:

UNIT HAS A DIAL CONTROL - SET TO MAX COLD. RECOMMEND DEFROST - ICE FORMING ON THE INNER CEILING OF THE UNIT.



Southeast Scientific
770.426.6678

Customer	CHOA - Center for Advanced Pediatrics		
Address	2174 N Druid Hills Road NE		
City, State, Zip	Atlanta, GA, 30329		
Ph	404-852-0765	Contact	Mimi Le
	1 st PM <input checked="" type="checkbox"/>		2 nd PM <input type="checkbox"/>

Unit Preventative Maintenance – Refrigerator / Freezer

Date of Service: 9/19/24 Next Service Due: 3/19/25 Room # 5126
 Model # HCUCFS-0404 Serial # VWR-051808050362-1809 Mfgr VWR
 SES # 9425 Customer # SHURI Tech CD
 NIST Meter Serial # T-281182 NIST Meter Calibration Due 6/26/26

	Pass	Fail	N/A	NEW		Yes	No	N/A
Door Seal	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Evap Drain Clear	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Hinges	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		CFM	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Battery	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Clean Condenser	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Gaskets	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Insulation Lines	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Handles	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					

Temperature

Setpoint D 2 Display NA NIST 3.6C Chart Recorder NA Adjust NA

Compressor

115 Voltage 0.6 Amperage

Inverter Normal
 No Inverter
 Inverter Error (See Notes)

1.1 Unit Total Amperage NA High Alarm Set NA Battery Charge Circuit
0.8 Amps Rated NA Low Alarm Set GOOD Wiring Connections
115 Volts Rated

Technician's Notes:

UNIT HAS A DIAL CONTROL
